

Home care for your implant retained auricular prosthesis

Step 1




Once a day, use the provided children's toothbrush and a mild liquid soap to clean the back side of the prosthesis. Use colorless hand soap with no strong fragrance, do not use antibacterial soap. Scrub the magnets or gold clips using the brush and soap lather.

Step 2



Using your fingers and a small amount of soap lather, wash the front side of the prosthesis. Rinse the prosthesis under cool or lukewarm water and gently dry with a cloth. Never use hot water or any extreme heat (i.e. sunlight) on the prosthesis.

Step 3



The front margins of the prosthesis are very thin (see photo). Careful attention must be made to not tear or damage this edge

Step 4



Store the prosthesis in the provided clean, dry container. Place the prosthesis in the container with the front side facing down. Doing this will keep the edges from curling

REC176-V02-10/25